



SLC23A1 Antibody, Biotin conjugated

Product Code	CSB-PA887038LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHI7
Immunogen	Recombinant Human Solute carrier family 23 member 1 protein (173-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 23 member 1 (Na(+)/L-ascorbic acid transporter 1) (Sodium-dependent vitamin C transporter 1) (hSVCT1) (Yolk sac permease-like molecule 3), SLC23A1, SVCT1 YSPL3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC23A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.