



SLC23A3 Antibody, FITC conjugated

Product Code	CSB-PA761358LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PIS1
Immunogen	Recombinant Human Solute carrier family 23 member 3 protein (391-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 23 member 3 (HPC E2-binding protein 3) (Na(+)/L-ascorbic acid transporter 3) (Sodium-dependent vitamin C transporter 3), SLC23A3, E2BP3 SVCT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC23A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.