



SLC25A26 Antibody, Biotin conjugated

Product Code	CSB-PA761690LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70HW3
Immunogen	Recombinant Human S-adenosylmethionine mitochondrial carrier protein (116-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-adenosylmethionine mitochondrial carrier protein (Mitochondrial S-adenosylmethionine transporter) (Solute carrier family 25 member 26), SLC25A26, SAMC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A26