



SLC25A4 Antibody, HRP conjugated

Product Code	CSB-PA021510LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12235
Immunogen	Recombinant Human ADP/ATP translocase 1 protein (40-74AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP/ATP translocase 1 (ADP,ATP carrier protein 1) (ADP,ATP carrier protein, heart/skeletal muscle isoform T1) (Adenine nucleotide translocator 1) (ANT 1) (Solute carrier family 25 member 4), SLC25A4, ANT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SLC25A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.