



SLC25A40 Antibody, FITC conjugated

Product Code	CSB-PA840561LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBP6
Immunogen	Recombinant Human Solute carrier family 25 member 40 protein (40-102AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 25 member 40 (Mitochondrial carrier family protein), SLC25A40, MCFP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.