




SLC25A46 Antibody

Product Code	CSB-PA856895LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AG3
Immunogen	Recombinant Human Solute carrier family 25 member 46 protein (1-102AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 25 member 46, SLC25A46, TB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A46

Image



Western blot

All lanes: SLC25A46 antibody at 6µg/ml

Lane 1: Mouse brain tissue


Lane 2: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 47, 31, 38 kDa

Observed band size: 47 kDa



Immunofluorescent analysis of HeLa cells using CSB-PA856895LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.