



SLC25A51 Antibody, Biotin conjugated

Product Code	CSB-PA875649LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1U9
Immunogen	Recombinant Human Solute carrier family 25 member 51 protein (1-35AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 25 member 51 (Mitochondrial carrier triple repeat protein 1), SLC25A51, MCART1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A51