




SLC35B4 Antibody

Product Code	CSB-PA853373LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969S0
Immunogen	Recombinant Human UDP-xylose and UDP-N-acetylglucosamine transporter protein (174-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-xylose and UDP-N-acetylglucosamine transporter (Solute carrier family 35 member B4) (YEA4 homolog), SLC35B4, YEA4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC35B4

Image



Immunofluorescent analysis of A2780 cells using CSB-PA853373LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.