



SLC35B4 Antibody, Biotin conjugated

Product Code	CSB-PA853373LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969S0
Immunogen	Recombinant Human UDP-xylose and UDP-N-acetylglucosamine transporter protein (174-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-xylose and UDP-N-acetylglucosamine transporter (Solute carrier family 35 member B4) (YEA4 homolog), SLC35B4, YEA4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC35B4