



SLC35B4 Antibody, Biotin conjugated

Product Code	CSB-PA853373LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969S0
Immunogen	Recombinant Human UDP-xylose and UDP-N-acetylglucosamine transporter protein (174-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-xylose and UDP-N-acetylglucosamine transporter (Solute carrier family 35 member B4) (YEA4 homolog), SLC35B4, YEA4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC35B4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.