



SLC35C1 Antibody, Biotin conjugated

Product Code	CSB-PA839285LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A29
Immunogen	Recombinant Human GDP-fucose transporter 1 protein (287-364AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDP-fucose transporter 1 (Solute carrier family 35 member C1), SLC35C1, FUCT1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SLC35C1