



# SLC35C1 Antibody, HRP conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA839285LB01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q96A29   |
| Immunogen           | Recombinant Human GDP-fucose transporter 1 protein (287-364AA)                   |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | HRP  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | GDP-fucose transporter 1 (Solute carrier family 35 member C1), SLC35C1, FUCT1    |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Tags & Cell Markers  |
| Target Names        | SLC35C1  |