



SLC35F5 Antibody, HRP conjugated

Product Code	CSB-PA840985LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV83
Immunogen	Recombinant Human Solute carrier family 35 member F5 protein (122-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 35 member F5 (Hepatitis C virus NS5A-transactivated protein 3) (HCV NS5A-transactivated protein 3), SLC35F5, NS5ATP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC35F5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.