



SLC35F6 Antibody, HRP conjugated

Product Code	CSB-PA843146LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N357
Immunogen	Recombinant Human Solute carrier family 35 member F6 protein (337-371AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 35 member F6 (ANT2-binding protein) (ANT2BP) (Transport and Golgi organization 9 homolog), SLC35F6, C2orf18
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC35F6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.