



SLC38A2 Antibody, HRP conjugated

Product Code	CSB-PA847210LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QD8
Immunogen	Recombinant Human Sodium-coupled neutral amino acid transporter 2 protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-coupled neutral amino acid transporter 2 (Amino acid transporter A2) (Protein 40-9-1) (Solute carrier family 38 member 2) (System A amino acid transporter 2) (System A transporter 1) (System N amino acid transporter 2), SLC38A2, ATA2 KIAA1382 SAT2 SNAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC38A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.