




SLC38A5 Antibody

Product Code	CSB-PA837853LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUX1
Immunogen	Recombinant Human Sodium-coupled neutral amino acid transporter 5 protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-coupled neutral amino acid transporter 5 (Solute carrier family 38 member 5) (System N transporter 2), SLC38A5, JM24 SN2 SNAT5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC38A5

Image



Immunofluorescent analysis of HeLa cells using CSB-PA837853LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.