



SLC38A5 Antibody, Biotin conjugated

Product Code	CSB-PA837853LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUX1
Immunogen	Recombinant Human Sodium-coupled neutral amino acid transporter 5 protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-coupled neutral amino acid transporter 5 (Solute carrier family 38 member 5) (System N transporter 2), SLC38A5, JM24 SN2 SNAT5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC38A5