



SLC38A7 Antibody, HRP conjugated

Product Code	CSB-PA889118HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVC3
Immunogen	Recombinant Human Putative sodium-coupled neutral amino acid transporter 7 protein (1-55AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative sodium-coupled neutral amino acid transporter 7 (Solute carrier family 38 member 7), SLC38A7, SNAT7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC38A7