



SLC38A7 Antibody, HRP conjugated

Product Code	CSB-PA889118HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVC3
Immunogen	Recombinant Human Putative sodium-coupled neutral amino acid transporter 7 protein (1-55AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative sodium-coupled neutral amino acid transporter 7 (Solute carrier family 38 member 7), SLC38A7, SNAT7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC38A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.