




SLC39A11 Antibody

Product Code	CSB-PA854023LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1S5
Immunogen	Recombinant Human Zinc transporter ZIP11 protein (93-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP11 (Solute carrier family 39 member 11) (Zrt- and Irt-like protein 11) (ZIP-11), SLC39A11, C17orf26 ZIP11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC39A11

Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA854023LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)