





SLC52A2 Antibody

Product Code	CSB-PA881008LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAB3
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 2 protein (217-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 2 (Porcine endogenous retrovirus A receptor 1) (PERV-A receptor 1) (Protein GPR172A) (Riboflavin transporter 3) (hRFT3), SLC52A2, GPR172A PAR1 RFT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC52A2

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA881008LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA881008LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.