



SLC52A2 Antibody, Biotin conjugated

Product Code	CSB-PA881008LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAB3
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 2 protein (217-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 2 (Porcine endogenous retrovirus A receptor 1) (PERV-A receptor 1) (Protein GPR172A) (Riboflavin transporter 3) (hRFT3), SLC52A2, GPR172A PAR1 RFT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC52A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.