



SLC6A16 Antibody, Biotin conjugated

Product Code	CSB-PA863929LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZN6
Immunogen	Recombinant Human Orphan sodium- and chloride-dependent neurotransmitter transporter NTT5 protein (1-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Orphan sodium- and chloride-dependent neurotransmitter transporter NTT5 (Solute carrier family 6 member 16), SLC6A16, NTT5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC6A16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.