




PRELID3B Antibody

Product Code	CSB-PA896708LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3B1
Immunogen	Recombinant Human PRELI domain containing protein 3B protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRELI domain containing protein 3B (Protein slowmo homolog 2), PRELID3B, C20orf45 SLMO2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRELID3B

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA896708LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.