



PRELID3B Antibody, Biotin conjugated

Product Code	CSB-PA896708LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3B1
Immunogen	Recombinant Human PRELI domain containing protein 3B protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRELI domain containing protein 3B (Protein slowmo homolog 2), PRELID3B, C20orf45 SLMO2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRELID3B