



# SMAD1 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA618998LC01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q15797  |
| <b>Immunogen</b>           | Recombinant Human Mothers against decapentaplegic homolog 1 protein (202-335AA)   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Mothers against decapentaplegic homolog 1 (MAD homolog 1) (Mothers against DPP homolog 1) (JV4-1) (Mad-related protein 1) (SMAD family member 1) (SMAD 1) (Smad1) (hSMAD1) (Transforming growth factor-beta-signaling protein 1) (BSP-1), SMAD1, BSP1 MADH1 MADR1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Signal Transduction   |
| <b>Target Names</b>        | SMAD1   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |