



SMIM14 Antibody, HRP conjugated

Product Code	CSB-PA857012LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QK8
Immunogen	Recombinant Human Small integral membrane protein 14 protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small integral membrane protein 14, SMIM14, C4orf34
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SMIM14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.