



SMPDL3A Antibody, HRP conjugated

Product Code	CSB-PA842119LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92484
Immunogen	Recombinant Human Acid sphingomyelinase-like phosphodiesterase 3a protein (51-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acid sphingomyelinase-like phosphodiesterase 3a (ASM-like phosphodiesterase 3a) (EC 3.1.4.-), SMPDL3A, ASML3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SMPDL3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.