



SMUG1 Antibody

Product Code	CSB-PA706636LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53HV7
Immunogen	Recombinant Human Single-strand selective monofunctional uracil DNA glycosylase protein (1-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-strand selective monofunctional uracil DNA glycosylase (EC 3.2.2.-), SMUG1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SMUG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.