



SNPH Antibody

Product Code	CSB-PA022317LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15079
Immunogen	Recombinant Human Syntaphilin protein (1-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaphilin, SNPH, KIAA0374
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNPH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.