



SNRNP25 Antibody, Biotin conjugated

Product Code	CSB-PA863129LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV90
Immunogen	Recombinant Human U11/U12 small nuclear ribonucleoprotein 25 kDa protein (1-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U11/U12 small nuclear ribonucleoprotein 25 kDa protein (U11/U12 snRNP 25 kDa protein) (U11/U12-25K) (Minus-99 protein), SNRNP25, C16orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRNP25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.