





SNX6 Antibody

Product Code	CSB-PA890766LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNH7
Immunogen	Recombinant Human Sorting nexin-6 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-6 (TRAF4-associated factor 2) [Cleaved into: Sorting nexin-6, N-terminally processed], SNX6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNX6

Image



Western blot
 All lanes: SNX6 antibody at 4µg/ml + HeLa whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 47, 34 kDa
 Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA890766LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.