



SNX6 Antibody, Biotin conjugated

Product Code	CSB-PA890766LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNH7
Immunogen	Recombinant Human Sorting nexin-6 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-6 (TRAF4-associated factor 2) [Cleaved into: Sorting nexin-6, N-terminally processed], SNX6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNX6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.