



Socs3 Antibody, Biotin conjugated

Product Code	CSB-PA022392LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O88583
Immunogen	Recombinant Rat Suppressor of cytokine signaling 3 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppressor of cytokine signaling 3 (SOCS-3) (Cytokine-inducible SH2 protein 3), Socs3, Cish3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Socs3