





SORBS3 Antibody

Product Code	CSB-PA022406LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60504
Immunogen	Recombinant Human Vinexin protein (1-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vinexin (SH3-containing adapter molecule 1) (SCAM-1) (Sorbin and SH3 domain-containing protein 3), SORBS3, SCAM1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SORBS3


Image



Western blot
 All lanes: SORBS3 antibody at 8µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 76, 37 kDa
 Observed band size: 76 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA022406LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA022406LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.