




SPANXC Antibody

Product Code	CSB-PA878921LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY87
Immunogen	Recombinant Human Sperm protein associated with the nucleus on the X chromosome C protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm protein associated with the nucleus on the X chromosome C (Cancer/testis antigen 11.3) (CT11.3) (Cancer/testis-associated protein CTP11) (Nuclear-associated protein SPAN-Xc) (SPANX-C) (SPANX family member C), SPANXC
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPANXC

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA878921LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.