



SPATS2L Antibody, Biotin conjugated

Product Code	CSB-PA868324HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUQ6
Immunogen	Recombinant Human SPATS2-like protein (2-488AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SPATS2-like protein (DNA polymerase-transactivated protein 6) (Stress granule and nucleolar protein) (SGNP), SPATS2L, DNAPTP6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPATS2L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.