



SPCS1 Antibody, HRP conjugated

Product Code	CSB-PA896543LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6A9
Immunogen	Recombinant Human Signal peptidase complex subunit 1 protein (70-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptidase complex subunit 1 (EC 3.4.-.-) (Microsomal signal peptidase 12 kDa subunit) (SPase 12 kDa subunit), SPCS1, SPC12
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPCS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.