



SPCS1 Antibody, HRP conjugated

Product Code	CSB-PA896543LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6A9
Immunogen	Recombinant Human Signal peptidase complex subunit 1 protein (70-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptidase complex subunit 1 (EC 3.4.-.-) (Microsomal signal peptidase 12 kDa subunit) (SPase 12 kDa subunit), SPCS1, SPC12
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPCS1