




SPERT Antibody

Product Code	CSB-PA843308LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NA61
Immunogen	Recombinant Human Spermatid-associated protein (1-412AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatid-associated protein (Protein chibby homolog 2), SPERT, CBY2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPERT

Image



Western blot
 All lanes: SPERT antibody at 7µg/ml
 Lane 1: Rat gonadal tissue
 Lane 2: Mouse lung tissue
 Lane 3: A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 52, 48 kDa
 Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.