




# SPERT Antibody

<b>Product Code</b>	CSB-PA843308LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8NA61
<b>Immunogen</b>	Recombinant Human Spermatid-associated protein (1-412AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Spermatid-associated protein (Protein chibby homolog 2), SPERT, CBY2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	SPERT

## Image



Western blot  
All lanes: SPERT antibody at 7µg/ml  
Lane 1: Rat gonadal tissue  
Lane 2: Mouse lung tissue  
Lane 3: A549 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 52, 48 kDa  
Observed band size: 52 kDa