




SPRTN Antibody

Product Code	CSB-PA863937LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H040
Immunogen	Recombinant Human SprT-like domain-containing protein Spartan protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SprT-like domain-containing protein Spartan (DNA damage protein targeting VCP) (DVC1) (Protein with SprT-like domain at the N terminus) (Spartan), SPRTN, C1orf124 DVC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SPRTN

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA863937LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.