



SRFBP1 Antibody, Biotin conjugated

Product Code	CSB-PA839864LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEF9
Immunogen	Recombinant Human Serum response factor-binding protein 1 protein (201-429AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum response factor-binding protein 1 (SRF-dependent transcription regulation-associated protein) (p49/STRAP), SRFBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SRFBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.