




ST6GALNAC6 Antibody

Product Code	CSB-PA850251LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969X2
Immunogen	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 6 protein (65-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 6 (EC 2.4.99.-) (GalNAc alpha-2,6-sialyltransferase VI) (ST6GalNAc VI) (ST6GalNAcVI) (hST6GalNAc VI) (Sialyltransferase 7F) (SIAT7-F), ST6GALNAC6, SIAT7F
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ST6GALNAC6

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA850251LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.