




STAC3 Antibody

Product Code	CSB-PA839373LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MF2
Immunogen	Recombinant Human SH3 and cysteine-rich domain-containing protein 3 protein (1-253AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 and cysteine-rich domain-containing protein 3, STAC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology

Image



Western blot
 All lanes: STAC3 antibody at 12µg/ml + Rat lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 42, 38, 21 kDa
 Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.