





STAR Antibody

Product Code	CSB-PA022798LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49675
Immunogen	Recombinant Human Steroidogenic acute regulatory protein, mitochondrial protein (64-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroidogenic acute regulatory protein, mitochondrial (StAR) (START domain-containing protein 1) (StARD1), STAR, STARD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAR

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA022798LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA022798LA01HU at dilution of 1:100