



STAR Antibody, FITC conjugated

Product Code	CSB-PA022798LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49675
Immunogen	Recombinant Human Steroidogenic acute regulatory protein, mitochondrial protein (64-285AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroidogenic acute regulatory protein, mitochondrial (StAR) (START domain-containing protein 1) (StARD1), STAR, STARD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.