




STK33 Antibody

Product Code	CSB-PA861179LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYT3
Immunogen	Recombinant Human Serine/threonine-protein kinase 33 protein (231-514AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 33 (EC 2.7.11.1), STK33
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STK33

Image



Western blot

All lanes: STK33 antibody at 4µg/ml

Lane 1: Mouse kidney tissue

Lane 2: Hela whole cell lysate


Lane 3: Jurkat whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58, 51 kDa

Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA861179LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA861179LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.