



STK38L Antibody, Biotin conjugated

Product Code	CSB-PA896491LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2H1
Immunogen	Recombinant Human Serine/threonine-protein kinase 38-like protein (212-464AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 38-like (EC 2.7.11.1) (NDR2 protein kinase) (Nuclear Dbf2-related kinase 2), STK38L, KIAA0965 NDR2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	STK38L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.