



STK38L Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA896491LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y2H1 |
| Immunogen | Recombinant Human Serine/threonine-protein kinase 38-like protein (212-464AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Serine/threonine-protein kinase 38-like (EC 2.7.11.1) (NDR2 protein kinase) (Nuclear Dbf2-related kinase 2), STK38L, KIAA0965 NDR2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | STK38L |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |